



### **GASM-001 Home Use Slimming Machine 4 Channels 8 Pads**

### **Contents:**

Manual, 8 Pads, 4 Elastic Belts, Gel, Bag.

#### Features:

- Adjustable time for the treatment
- Selectable range of frequency
- Continuous & intermittent mode
- Relax & contraction time setting
- External terminal to operate on alternate supply of 12 volts dc.
- No strenuous exercises, No medicine, No crash dieting, No side effects.



- Output: 4 Channels 8 Pads
- Modes: Weight loss, Inch loss
- Digital Timer: 1-60 Minutes
- Power supply from 90 volts to 230 volts AC 50 Hz









### GASM-002 Home Use Slimming Machine 6 Channels 12 Pads

### Contents:

Includes Manual, 12 Pads, 6 Elastic Belts, Gel, Bag

### **Features:**

- Adjustable time for the treatment
- Selectable range of frequency
- Continuous & intermittent mode
- Relax & contraction time setting
- External terminal to operate on alternate supply of 12 volts dc.
- No strenuous exercises, No medicine, No crash dieting, No side effects

### Specifications:

- Output: 6 Channels 12 Pads
- Modes: Weight loss, Inch loss, Skin tightening
- Digital Timer:1-60 Minutes
- Power supply from 90 volts to 230 volts AC 50 Hz





















### Contents:

Includes Manual, 16 Pads, 8 Elastic Belts, Gel, Bag

### - Adimetalala tima familia a turatura ut

- Adjustable time for the treatment
- Selectable range of frequencyContinuous & intermittent mode
- Relax & contraction time setting
- neida di contraction time setti
- Made of Steel Body
- External terminal to operate on alternate supply of 12 volts dc.
- No strenuous exercises, No medicine, No crash dieting, No side effects

### **Specifications:**

- Output: 8 Channels 16 Pads
- Modes: Weight loss, Inch loss, Skin tightening
- Digital Timer:1-60 Minutes
- Power supply from 90 volts to 230 volts AC 50 Hz







### **Contents:**

Includes Manual, 16 Pads, 8 Elastic Belts, Gel, Bag

#### **Features:**

- Adjustable time for the treatment
- Selectable range of frequency
- Continuous & intermittent mode
- Relax & contraction time setting
- Made of Acrylic Housing
- External terminal to operate on alternate supply of 12 volts dc.
- No strenuous exercises, No medicine, No crash dieting, No side effects

### **Specifications:**

- Output: 8 Channels 16 Pads
- Modes: Weight loss, Inch loss, Lifting, Figure correction, Toning, Firming, Massaging, Body shaping
- **Digital Timer:**1-60 Minutes
- Power supply from 90 volts to 230 volts AC 50 Hz









### **GASM-005 Body Shaping Machine 12 Channels 24 Pads**

### **Contents:**

Includes Manual, 24 Pads, 12 Elastic Belt, 1 Power Card, Bag

#### **Features:**

- · Adjustable time for the treatment
- Selectable range of frequency
- Continuous & intermittent mode
- Relax & contraction time setting
- External terminal to operate on alternate supply of 12 volts dc.
- No strenuous exercises, No medicine, No crash dieting, No side effects

### **Specifications:**

- Output: 12 Channels 24 Pads
- Modes: Weight loss, Inch loss, Toning
- Digital Timer:1-60 Minutes
- Power supply from 90 volts to 230 volts AC 50 Hz







## GASM-006 Full Body Slimming Machine 12 Channels 24 Pads LED Programs

- Adjustable time for the treatment
- Selectable range of frequency
- Continuous & intermittent mode
- Relax & contraction time setting
- External terminal to operate on alternate supply of 12
- No strenuous exercises, No medicine, No crash dieting, No side effects

### **Specifications:**

Features:

- Output: 12 Channels 24 Pads LED Programs
- Modes: Weight loss, Inch loss, Lifting, Figure correction, Toning, Firming, all in one, Body shaping
- **Digital Timer:**1-60 Minutes
- Power supply from 90 volts to 230 volts AC 50 Hz

































### **Features:**

- Adjustable time for the treatment
- Selectable range of frequency
- Continuous & intermittent mode
- Relax & contraction time setting
- External terminal to operate on alternate supply of 12 volts dc.
- No strenuous exercises, No medicine, No crash dieting, No side effects

### **Specifications:**

- Output: 12 Channels 24 Pads, 3 Programs
- Modes: Weight loss, Inch loss, Lifting, Figure correction, Toning, Firming, all in one, Body shaping
- Digital Timer: 1-60 Minutes
- Power supply from 90 volts to 230 volts AC 50 Hz

### GASM-008 Body Shaping Machine 16 Channels 32 Pads LED Machines

#### **Features:**

- Adjustable time for the treatment
- Selectable range of frequency
- Continuous & intermittent mode
- Relax & contraction time setting
- External terminal to operate on alternate supply of 12 volts dc.
- No strenuous exercises, No medicine, No crash dieting, No side effects

### **Specifications:**

- Output: 16 Channels 34 Pads
- Modes: Weight loss, Inch loss, Lifting, Figure correction, Toning, Firming, Body shaping
- **Digital Timer:**1-60 Minutes
- Power supply from 90 volts to 230 volts AC 50 Hz














### GASM-009 Body Shaping Machine 20 Channels 40 Pads LED Machines



### **Features:**

- Adjustable time for the treatment
- Selectable range of frequency
- Continuous & intermittent mode
- Relax & contraction time setting
- External terminal to operate on alternate supply of 12 volts dc.
- No strenuous exercises, No medicine, No crash dieting, No side effects

### **Specifications:**

- Output: 20 Channels 40 Pads
- Modes: Weight loss, Inch loss, Lifting, Figure correction, Toning, Firming, Body shaping
- Digital Timer: 1-60 Minutes
- Power supply from 90 volts to 230 volts AC 50 Hz













### **Steam Bath Machine**

### **GASM-0010 Steam Bath Machine**

Our company specialize in offering qualitative range of Steam Bath Equipment. These products are being offered to our customers at cost effective prices. Our vendors manufacture these products using high grade raw material and latest technology & machinery. The range of our products helps in relaxing the pores and tired muscles. Moreover, these products effectively cure acne, eczema, psoriasis, cellulite and obesity.

#### **Features:**

- · Light weight
- Folding waterproof cabin
- With remote control

### **Specifications:**

Voltage : 110-220 vPower : 77 w

## **Heating Therapy Slimming Machine**

### **GASM-0011 Cellotherm Heat Machine**

### **Features:**

- Adjustable time for the treatment
- · Selectable range of frequency
- Continuous & intermittent mode
- Relax & contraction time setting
- External terminal to operate on alternate supply of 12 volts dc.
- No strenuous exercises, No medicine, No crash dieting, No side effects

### **Specifications:**

- Output: 5 Pads, I Tummy, 2 Leg, 2 Arms
- Digital Timer:1-60 Minutes
- Power supply from 90 volts to 230 volts AC 50 Hz.



**Slimming Massager** 

### **GASM-0012 Thrive Tummy Tuck G 5 Massager**

#### Features:

- G5 massage assists with the appearance of lumpy and congested skin.
- G5 works on breaking down fatty deposits, cellulite and dispersing this with excess fluid through the lymphatic areas.
- This process allows an instant inch loss and a detox effect, as well as improvement in skin elasticity, collagen flow and blood circulation.
- G5 provide toning and firming of the skin, a reduction in cellulite and fatty tissues all resulting in temporary but measurable reduction in body measurements making those "too tight" clothes more comfortable.
- A wonderful improvement in cellulite and texture of skin and elasticity and detox result will be felt and seen with instant inch loss.

### 1 sitting reduces .5 to 1.5 inches

• No strenuous exercises, No medicine, No crash dieting, No side effects.
















 G10 massage assists with the appearance of lumpy and congested skin.

 G10 works on breaking down fatty deposits, cellulite and dispersing this with excess fluid through the lymphatic areas.


 This process allows an instant inch loss and a detox effect, as well as improvement in skin elasticity, collagen flow and blood circulation.

 G10 provide toning and firming of the skin, a reduction in cellulite and fatty tissues all resulting in temporary but measurable reduction in body measurements making those "too tight" clothes more comfortable.

 A wonderful improvement in cellulite and texture of skin and elasticity and detox result will be felt and seen with instant inch loss.

### 1 sitting reduces .5 to 1.5 inches

 No strenuous exercises, No medicine, No crash dieting, No side effects



- · Detoxifies and burns fat.
- Quickens lipolysis. Remove unwanted fats in your Tummy, waist, hip and four limbs
- Extracts waste from the body such as heavy metals, poisons etc.; activates hormones
- Speeds up metabolism and enhances absorption of nutrients.
- Acts as physical therapy for myalgia & neurosis.
  Balances the autonomic nervous system and has a beneficial effect on sore muscles, lactic acid and the
- Improves the immune system and prevents growth of cancer cells.
- Diminishes inflammation.
- Speeds up metabolism so that improves functioning of organs, such as liver, intestine and stomach
- Eliminate intense and pressure of life and promote working efficiency
- Recuperates insufficiency physiological exercise, promotes sleep and sleepiness

## GASM-0015 Foot Massager

#### Features:

- Massages entire lower leg area, including angles and feet
- Customizable massage with dual-mode kneading and vibration programs
- Three speed levels and three intensity settings
- 15-minute auto shut-off timer
- Adjustable tilt angle for easy and comfortable seating position
- Unit rotates to become a comfortable footrest
- Includes two removable fabric covers for easy cleaning





**Heat Blanket** 



















### **Cryopolysis Machine**

### GASM-0016 Cryolipolysis 3 in 1 Program/ GASM-0017 Cryopolysis

### Features:

Instant Fat Loss and inch-loss by Freezing and suction Technology with advance RF and Cavitation.

#### Cryolipolysis Treatment.

- Usually 30-60 minutes for each part (depending on thickness of fat and Health status, After 20-30 minutes of treatment, press the pause button to suspend the suction, check and add antifreeze fluid..)
- Can only treat Cellulite part by part, for example, to eliminate Cellulite of the entire waist, you need to complete several Courses.
- Take a another cooling procedure at the same position 20 to 25 days after the first course, which can reduce the extra fat layer thickness, twice treatment, can enhance the effect of cryolipolysis, but not recommended to accept treatment repeatedly at random.
- Cryolipolysis cool-sculpturing treatment is designed for those physically healthy people who have local fat deposits, but not suitable for the one with overall obesity. It can target at te fat at a specific part, such as waist, abdomen, back etc. The treatment is not a panacea, it cannot treat obesity. It is not a substitute for liposuction or other invasive treatments. Substantial weight loss should be achieved through a balanced diet and regular exercise, cryolipolysis cool-sculpturing treatment can only eliminate the extra fat.
- Use appropriate and constant suction intensity, the suction intensity is adjustable according to the patient's individual vulnerability. At the same time, the mode of suction and massage can be used.
- Qualified antifreeze fluid is necessary to avoid forstbite during the treatment process when the temperature is reduced to 0 to -40 degrees.
   (Consumables, Purchased by the user)
- After the treatment, skin injury, minor frostbite, or pigment may occur, pay attention to take good care of it.

### **Features**

The face, neck anti-aging. Shrink pores, whitening. The reduction the orange peel organization. Tightening skin, strengthening. Promotion metabolism, blood circulation. Replenishing moisture content. Decrease the excessive water and piled fat. Removing Stomach wrinkle. Complementing collagen, recovering resilience.